



In a demanding and rapidly changing healthcare environment, effective communications is the lifeline for doctors and nurses. The need to deliver efficient quality care is vital—hence the requirement for the latest, most reliable technology. At Telecor, we seek to address the communication needs of today's healthcare facilities. Our Nurse Courier PLUS Software offers high-quality, versatile communications, featuring Wireless Telephony and Pocket Pager integration, Wireless Device Assignments, Management Software and so much more.

Add this to the cutting-edge equipment of the Nurse Courier PLUS System—including Consoles, Nurse Call stations, Public Address and PBX integration—and you've got the ideal communications solution for the Healthcare Market.



NURSE COURIER® PLUS SOFTWARE


The Heart of Healthcare Performance



NURSE COURIER PLUS SOFTWARE

A PRACTICAL SOLUTION FOR HEALTHCARE FACILITIES

Nurse Courier PLUS Software is a modular application that integrates the Nurse Courier 3 Communication System to wireless technology to provide instant call notification using text messaging. With the benefit of wireless mobility, doctors, nurses and other support staff have access to real-time information anywhere in the facility. The result? Superior response times with improved patient care.



WIRELESS COMMUNICATION

- > Seamless integration with virtually any Wireless Phone or Pocket Pager Network
- > Routing of patient calls directly to wireless telephones and pocket pagers, promoting staff mobility
- > Customized tag messages for wireless devices, informing staff of specific patient requests
- > Routing of calls to personal cell phones, or email messaging to locations outside the facility, ensuring key personnel are never out of touch

EASY-TO-USE GRAPHICAL USER INTERFACE

- > Incorporate a graphic to represent the physical location of different activities
- > Assign and reassign call points to wireless devices with a simple "drag and drop" action
- > Forward calls to other wireless devices when staff are busy or on break
- > Program assignment schedules to automatically be reassigned to accommodate staffing requirements for different times of the day
- > Customize a series of escalation responses for unanswered calls
- > Send instant or scheduled text messages to wireless devices
- > Monitor call-point activity in a particular area on an "active alarms" screen
- > Activate "virtual call-points" through wireless phones or PC's, which can notify other wireless devices of events or emergencies

RESOURCEFUL MANAGEMENT SOFTWARE

- > View, on a PC, central reporting of staff and resident/patient activity in real-time
- > Generate comprehensive activity reports, such as total calls placed, average call response time, average presence in room time, total service requirements, and average staff response time
- > Instantly display call and patient information on Electronic Wallboards
- > Monitor multiple independent systems on a single point of control
- > Tailor the Nurse Courier 3 System and any wireless network configuration to fit a facility's exact needs

www.telecor.com

In the USA — 165 Kirkland Circle, Oswego, Illinois 60543 • Phone: (630) 236-8400 • Fax: (630) 978-2201
 In Canada — 1114 Westport Crescent, Mississauga, Ontario L5T 1G1 • Phone: (905) 564-0801 • Fax: (905) 564-0806
 In the UK — Unit 3a, 21 Coopers Court, Coopers Yard, Newport Pagnell; Milton Keynes MK16 8JS; Buckinghamshire • Phone: 44-(0)-1908-211782 • Fax: 44-(0)-1908-216946

All product information subject to change without notice.
 Copyright © 2004 by Telecor Incorporated.
 Telecor, the Telecor logo, and Nurse Courier are Registered Trademarks of Telecor Inc.

Rev: 1.0
 ETS- 2367
 Printed in Canada